

## DAX Copilot commands cheat sheet

## Useful commands for DAX Copilot on the mobile application

What you can say	What it does
Next field	Jump to the next note section
Previous field	Jump to the previous note section
Go to [SECTION]	Jump to the stated section within the note
Go to social history, Go to allergies, Go to review of systems	
Select	Select the word spoken
Select [xxxx] through [xxxx]	Select the specified text
Delete that	Delete the selected text
New paragraph	Begin a new paragraph

## Useful commands for DAX Copilot on the desktop

What you can say	What it does
Open DAX Copilot Show DAX Copilot Launch DAX Copilot	Launch the DAX Copilot window
Close DAX Copilot	Closes the DAX Copilot window
Provide DAX Copilot feedback Send DAX Copilot feedback Provide DAX feedback Send DAX feedback	Opens DAX Copilot in-app feedback
Transfer [SECTION] (e.g., "transfer history of present illness")	Transfer the content of the stated section into your EHR
Transfer selected	Transfer the content of the selected sections into the EHR
Transfer all	Transfer the content of all non-blank sections into EHR.
Delete that	Delete the selected text



## Useful Dragon Medical One commands

DAX Copilot also supports most general Dragon Medical One commands:

- Scratch that/Delete that (Deletes the last dictated text)
- Delete <XYZ> (XYZ is the word you want to delete)
- Select <XYZ> (XYZ is the word you want to select)
- Select that (selects last utterance)
- New line / new paragraph
- Insert before/after <XYZ>
- Select all
- Select <word> through <word>
- Unselect that
- Cap that (first letter in each word is capitalized)
- All caps that (all letters in selection are capitalized)

**Note:** To see a list of all available commands, say What Can I Say.

MEDICAL DEVICE DISCLAIMER. Nuance Products are not designed, intended, or made available as a medical device(s), and are not designed or intended to be a substitute for professional medical advice, diagnosis, treatment, or judgment and should not be used to replace or as a substitute for professional medical advice, diagnosis, treatment, or judgment. Nuance does not warrant that Nuance products will be sufficient for any medical purposes or meet the health or medical requirements of any person.